

Pluronic® F-127

P-3000 Pluronic® F-127 *20% solution in DMSO*
P-6866 Pluronic® F-127 *10% solution in water* *0.2 µm filtered*
P-6867 Pluronic® F-127 *low UV absorbance*

Introduction

Pluronic® F-127 is a nonionic, surfactant polyol (molecular weight approximately 12,500 daltons) that has been found to facilitate the solubilization of water-insoluble dyes and other materials in physiological media.¹ Pluronic F-127 has been used to help disperse acetoxymethyl (AM) esters of fluorescent ion indicators such as fura-2, indo-1, fluo-3 and SBFI; it appears to be required with SBFI-AM or PBFI-AM, and optional with other indicators. Pluronic F-127 may also be useful for dispersing other lipophilic probes. Appropriate controls should be performed to make certain that Pluronic F-127 is not altering the membrane properties of the cell. For the convenience of our customers, Molecular Probes offers Pluronic F-127 in three forms: 1 mL of a 20% (w/v) solution in DMSO (P-3000), 30 mL of a 0.2 µm–filtered 10% (w/v) solution in water (P-6866) and 2 g of the solid (P-6867).

Storage and Handling

Dissolve 2 g of Pluronic F-127 in 10 mL of anhydrous dimethylsulfoxide (DMSO) to give a 20% (w/v) stock solution. This may require heating at ~40° for about 20 minutes. Store both solid and solution at room temperature. Do not refrigerate or freeze the solutions because Pluronic F-127 may come out of solution. If the product does crystallize, it can be reheated until it goes back into solution.

Application

The experimental conditions for loading cells with AM esters varies with cell type due to differences both in probe uptake and in the intracellular esterase activity required for hydrolysis of the AM esters. Solutions of the AM esters in DMSO must be kept anhydrous since the solvent will readily take up moisture, leading to loss of cell-loading efficacy. Pluronic F-127 should be added only to working solutions. Typically, a small volume of the AM ester, dissolved at 1–5 mM in DMSO, is mixed with the 20% (w/v) Pluronic F-127 stock solution in DMSO at a ratio of 1:1 immediately before use. The solution of AM ester and Pluronic F-127 is then diluted into the cell-loading buffer to achieve a final AM ester concentration of between 1 µM and 10 µM and the cells are incubated for between 10 minutes and 1 hour. The final concentration of Pluronic F-127 is normally kept at or below 0.1%. More weakly fluorescent indicators, such as the AM esters of SBFI, PBFI, quin-2 and Fura Red may require more concentrated loading solutions and correspondingly greater amounts of Pluronic F-127.² In general it is desirable to use the minimum amount of AM ester needed to achieve adequate fluorescence signal to noise levels. Loading may be done at any temperature that is tolerable for the cells. Note that the incubation temperature generally affects the extent of intracellular dye compartmentalization.^{3,4} After labeling, the cells are washed with fresh medium before beginning the experiment.

References

1. J Membrane Biol 19, 1 (1974); 2. J Biol Chem 265, 19543 (1990); 3. Methods Enzymol 302, 341 (1999); 4. Methods Enzymol 307, 441 (1999).

Product List *Current prices may be obtained from our Web site or from our Customer Service Department.*

Cat #	Product Name	Unit Size
P-6867	Pluronic® F-127 *low UV absorbance*	2 g
P-3000	Pluronic® F-127 *20% solution in DMSO*	1 mL
P-6866	Pluronic® F-127 *10% solution in water* *0.2 µm filtered*	30 mL

Contact Information

Further information on Molecular Probes' products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site — www.probes.com — for the most up-to-date information

Molecular Probes, Inc.

PO Box 22010, Eugene, OR 97402-0469

Phone: (541) 465-8300 • Fax: (541) 344-6504

Customer Service: 7:00 am to 5:00 pm (Pacific Time)

Phone: (541) 465-8338 • Fax: (541) 344-6504 • order@probes.com

Toll-Free Ordering for USA and Canada:

Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Assistance: 8:00 am to 4:00 pm (Pacific Time)

Phone: (541) 465-8353 • Fax: (541) 465-4593 • tech@probes.com

Molecular Probes Europe BV

PoortGebouw, Rijnsburgerweg 10

2333 AA Leiden, The Netherlands

Phone: +31-71-5233378 • Fax: +31-71-5233419

Customer Service: 9:00 to 16:30 (Central European Time)

Phone: +31-71-5236850 • Fax: +31-71-5233419

eurorder@probes.nl

Technical Assistance: 9:00 to 16:30 (Central European Time)

Phone: +31-71-5233431 • Fax: +31-71-5241883

eurotech@probes.nl

Molecular Probes' products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Several of Molecular Probes' products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to busdev@probes.com. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2001, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.